

Chest Freezer

- Temperature Range -10°C to -25°C, -10°C to -40°C.
- Digital temperature controller with display.
- Pt 1000 sensor with 0.1° accuracy on the display.
- Painted steel external and aluminium internal.
- 40mm foaming insulation.
- 1 unit steel wire basket is included.
- Delay start function to protect compressor.
- Rapid cool down to -25°C in 2-3 hours
- Mounted on castors.
- Door with safety locker.
- CFC Free Refrigerant R134a
- Certification: ISO 9001, ISO 14001, ISO 13485
- Power supply: 220V 50Hz (110V 60Hz is available)
- Warranty period: 2 years on all parts.



ST-CF40-251

Specification of Chest Freezer

Model	Capacity (L)	External (W x D x H) mm	Internal (W x D x H) mm	Power Rating (W)	Weight (kg)
ST-CF25-110	110	549 x 549 x 845	410 x 410 x 654	145	46
ST-CF25-166	166	556 x 906 x 937	430 x 780 x 480	160	71
ST-CF25-196	196	556 x 1056 x 937	430 x 930 x 480	180	67
ST-CF25-226	226	556 x 1206 x 937	430 x 1080 x 480	207	71
ST-CF25-358	358	730 x 1204 x 968	530 x 1080 x 625	320	97
ST-CF25-508	508	730 x 1554 x 968	530 x 1400 x 685	375	115
ST-CF40-110	110	583 x 619 x 870	410 x 410 x 654	260	58
ST-CF40-251SS	251	730 x 1204 x 968	430 x 1050 x 556	265	90
ST-CF40-351	351	730 x 1554 x 968	430 x 1300 x 627	300	110

NOTE: -10°C to -40°C featured stainless steel interior for 251L and 351L models. All equipped with alarm system for high/low temperature against set temperature.