

## Digital Constant Temperature Hotplate



- Non-warping steel hotplate with heating elements embedded throughout the surface produces even heat distribution.
- Stainless steel 304 outer casing.
- PID digital temperature controller displays set and actual temperature, maximum working duration: 48 hours.
- Working temperature up to 350°C.

### Specification of Digital Constant Temperature Hotplate

Model	ST-HP120S	ST-HP150S	ST-HP200S	ST-HP300S
Temperature Range	120 x 120	150 x 150	200 x 300	300 x 400
Plate Size (mm)	Up to 350°C			
Power Supply	AC 220V 50Hz			

## Circular Hotplate



- Electronic controller.
- Heating elements embedded throughout the non-warping steel surface produces even heat distribution.
- Power-on indicator.
- Working temperature up to 400°C.

### Specification of Circular Hotplate

Model	ST-HP125R	ST-HP185R
Plate Size (mm)	Ø 125mm	Ø 185mm
Temperature Range	Up to 400°C	
Power Supply	AC 220V 50Hz	