



PROPISEL RACKS

Made of high density white polypropylene. Easy to assemble. Stackable. Autoclavable. These racks do not float and are suitable for water baths for temperatures up to 140 °C. Each rack has a cross reference numbering systems for easy tube identification. Easily cleaned with a disinfectant. Colour white. Accessory rack identification cards are available in: yellow, blue and red.

Models

Part No.	Number of tubes	For tubes up to Ø mm	External dims. cm		
			Height	Width	Depth
1001115	90	13	7	10.5	25
1001116	60	16	7	10.5	25
1001117	40	20	7	10.5	25
1001118	24	25	7	10.5	25
1001119	21	30	7	10.5	25

ACCESSORY

Identification cards, autoclavable.

Part No.

1001120	Yellow
1001121	Blue
1001122	Red



PROPISEL-FIX RACKS

Made of high density polypropylene. Comes with decanting aid that holds the tubes in place in the rack while the complete rack is emptied in one operation. Features: Same feature as the "Propisel".

Image showing the decanting aid in use and tubes still in place.

Models


Part No.	Number of tubes	For tubes up to Ø mm	External dims. cm		
			Height	Width	Depth
1001123	90	11 to 13	7	10.5	25
1001124	60	14 to 16	7	10.5	25
1001125	40	19 to 20	7	10.5	25

ACCESSORY

Identification cards, autoclavable.

Part No.

1001120	Yellow
1001121	Blue
1001122	Red



RACK FOR 20 10 mm □ SQUARE CUVETTES

Made of polypropylene. Positioning cross-reference. Height 35 mm. Width 70 mm. Depth 210mm. Colour: white. Part No. 1001528



CIRCULAR RACK

Made in polyethylene for 4 tubes of 25 mm Ø, 8 tubes of up to 19 mm Ø and 16 up side down peg tube supports. Dimensions: 95 high x 180 mm Ø. Part No. 1001339



RACK FOR 1.5 AND 2 ml MICROTUBES

Made of Polycarbonate, stackable, numbered positions and autoclavable. Number of tubes: 24 microtubes. Diameter: up to 11 mm Ø. Part No. 1001340



RACK FOR 1.5 AND 2 ml MICROTUBES

Made of solid polypropylene, ideal for refrigerator and freezer storage. Numbered positions and autoclavable. Capacity: up to 80 micro-tubes. Part No. 1001341



RACK FOR 20 MICROTUBES OF 1.5 TO 2 ml

Made of polypropylene. Racks can be linked and fitted together. Numbered locations. Total length 21 cm. Part No. 1001523