

Low Speed Self-balance Centrifuge

- Microprocessor controller, LCD display.
- Maintenance free DC brushless motor.
- Stainless steel chamber with electronic interlocking lid.
- Fast acceleration and deceleration less than 30 seconds for acceleration and 35 seconds in deceleration.
- Storage memory for 10 programmes.
- Short spin, timer 0-99 minutes, special button for Rcf display.



Specification of Low Speed Self-balance Centrifuge

Model		ST-CT500
Maximum Speed		5000 r.p.m.
Maximum R.c.f.		4100xg
External Dimension (L x W x H) mm		370 x 432 x 250mm
Weight		15kg
Accessories	ST-CT500-1	Angle rotor 24 x 10ml, 4000r.p.m., 2580xg
	ST-CT500-2	Swing out rotor 12 x 15ml, 5000r.p.m., 4100xg
	ST-CT500-3	Swing out rotor 16 x 15ml, 5000r.p.m., 4100xg
	ST-CT500-4	Swing out rotor 24 x 5ml, 5000r.p.m., 4100xg

High Speed Refrigerated Centrifuge

- Microprocessor controller, LCD display.
- Maintenance free DC brushless motor.
- Electronic door interlock for safety operation.
- Special button for R.c.f. display.
- Inclusive of rotor 12 x 1.5/ 2.2ml.



Specification of High Speed Refrigerated Centrifuge

Model		ST-CT160
Maximum Speed		16000 r.p.m.
Maximum R.c.f.		17800xg
Temperature Setting Range		-20°C ~ 40°C
Refrigeration System		R404a, non-freon refrigeration system
Timer		0-99 min
External Dimension (L x W x H) mm		400 x 560 x 300
Weight		65kg
Power Supply		AC 220V 50Hz
Accessories	ST-CT160-1	Angle rotor 12 x 1.5/ 2.2ml, 16000rpm, 17800xg
	ST-CT160-2	Angle rotor 10 x 5ml, 13000rpm, 11400xg
	ST-CT160-3	Angle rotor 18 x 0.5ml, 16000rpm, 14800xg
	ST-CT160-4	Angle rotor 24 x 1.5/ 2.2ml, 13000rpm, 15100xg