


Visible Spectrophotometer



| Visible Spectrophotometer | |
|---------------------------|--|
| Model | ST-722 |
| Optical System | Single beam, 1200 lines/ mm |
| Wavelength Range | 325-1000nm |
| Spectra Bandwidth | 6nm |
| Wavelength Accuracy | ± 2nm |
| Reproducibility | ≤1nm |
| Wavelength Readout | 4 digit LED display |
| Stray Light | ≤0.5% (T) at 360nm |
| Photometric Range | 0-125%T, -0.097-1.99A, 0-1999C (0-1999F) |
| Photometric Accuracy | ± 1%T |
| Light Source | Tungsten-halogen lamp |
| Detector | Photocell |
| Data Output | RS 232 serial port |
| Printer Output | Serial port |
| Voltage | 85-260V, 50Hz |
| Dimension (L x W x D) mm | 430 x 310 x 200 |
| Working Environment | 5~35°C, below 80% RH |
| Net Weight | 8kg |
| Standard Accessories | 4 glass cells, dust cover and operation manual |

UV/ Vis Spectrophotometer



| UV/ Vis Spectrophotometer | |
|---------------------------|--|
| Model | ST-7504C |
| Optical System | Single beam, CT monochromater |
| Wavelength Range | 200-1000nm |
| Spectra Bandwidth | 4nm |
| Wavelength Accuracy | ± 2nm |
| Wavelength Repeatability | ≤1.0nm |
| Stray Light | < 3% T |
| Photometric Range | 0-125%T, -0.097-2.50A, 0-1999C (0-1999F) |
| Photometric Accuracy | ± 0.5%T |
| Photometric Repeatability | < 0.2%T |
| Light Source | Tungsten-halogen lamp & deuterium |
| Display | LCD, 2 x 20 |
| Detector | Silicon Photocell |
| Data Output | RS 232C |
| Dimension (L x W x D) mm | 440 x 370 x 200 |
| Working Environment | 5~35°C, below 80% RH |
| Net Weight | 17kg |
| Standard Accessories | 4 glass cuvette 10mm, 2 quartz cuvette 10mm & dust cover |

Flake Ice Maker

- Compact design, attractive appearance.
- Fully automatic system air cool compressor system ventilate from front to rear.
- Safety system for overloading power and water shortage.
- High efficiency CFC free compressor, economical and low energy consumption.



| Specification of Flake Ice Maker | | | |
|-----------------------------------|-----------------------|-----------------|-----------------|
| Model | ST-FIM50 | ST-FIM100 | ST-FIM200 |
| Capacity | 50kg | 100kg | 200kg |
| Storage Bin | 10kg | 20kg | 35kg |
| Cabinet | Stainless steel | | |
| External Dimension (W x D x H) mm | 398 x 542 x 732 | 548 x 612 x 872 | 635 x 612 x 945 |
| Packing Dimension (W x D x H) mm | 465 x 615 x 790 | 615 x 685 x 932 | 712 x 690 x 994 |
| Refrigerant | R134a | | |
| Insulation Material | CFC free polyurethane | | |
| Net Weight | 40kg | 58kg | 90kg |