

Floor Model & Table Top Drying Oven

- Interior of oven is made of stainless steel, and exterior is made of high-quality steel.
- Microprocessor PID temperature controller with temperature control protection, digital display and timing function ensures a precise reliable temperature control.
- Hot air circulating system and air ducts to ensure a uniform distribution and high temperature in the chamber, back horizontal air flow.
- Independent alarm system for over temperature.



Specification of Drying Oven

Model	Ambient +10 to 250°C	ST-70A	ST-140A	ST-240A	ST-420A	ST-620A	ST-920A
	Ambient +10 to 300°C	ST-70B	ST-140B	ST-240B	ST-420B	ST-620B	ST-920B
Capacity		70L	136L	220L	420L	620L	1000L
Interior Dimension (W x D x H) mm		400 x 350 x 500	450 x 550 x 550	500 x 600 x 750	600 x 550 x 1300	800 x 600 x 1300	1000 x 600 x 1600
Exterior Dimension (W x D x H) mm		600 x 550 x 755	600 x 710 x 845	880 x 720 x 930	745 x 750 x 1670	945 x 800 x 1690	1140 x 800 x 2150
Shelves (Adjustable)		2 (PCS)	2 (PCS)	2 (PCS)	3 (PCS)	4 (PCS)	5 (PCS)
Material	Inner: Stainless Steel; Outer: High Quality Steel, Thickness 1-1.5mm						
Door	Double safety glass window door						
Controller	Microprocessor PID multi-function controller						
Display	LED 4 digits display						
Timer	0-9999 minutes						
Circulation	Air inlet can be switch, low noise blower fan						
Temperature Precision	0.1°C						
Temperature Stability	± 1.0°C						
Input Power	ST-.....A	1050W	1500W	2100W	4000W	4500W	6000W
	ST-.....B	1240W	1740W	2500W	4500W	4800W	7500W
Net Weight		54 kg	64 kg	92 kg	150 kg	172 kg	196 kg
Electrical Requirement	AC 220V 50/60Hz						

Tabletop Drying Oven

- Interior of oven is made of stainless steel, and exterior is made of high-quality steel.
- Microprocessor PID temperature controller with temperature control protection, digital display and timing function ensures a precise reliable temperature control.
- Hot air circulating system and air ducts to ensure a uniform distribution and high temperature in the chamber, back horizontal air flow.
- Independent alarm system for over temperature.



Specification of Tabletop Drying Oven

Model	ST-23A	ST-53A	ST-123A	ST-203A
Capacity	23L	50L	105L	198L
Interior Dimension (W x D x H) mm	300 x 300 x 270	420 x 350 x 350	550 X 350 X 550	600 x 550 x 600
Exterior Dimension (W x D x H) mm	500 x 450 x 440	710 x 520 x 520	840 x 530 x 730	880 x 765 x 780
Shelves	Adjustable 2 or 3 shelves			
Material	Inner: Stainless Steel; Outer: High Quality Steel, Thickness 1-1.5mm			
Door	Double safety glass window door			
Controller	Microprocessor PID multi-function controller			
Display	LED 4 digits display			
Timer	0-9999 minutes			
Circulation	Air inlet can be switch, low noise blower fan			
Temperature Range	Ambient +10 to 200°C			
Temperature Precision	0.1°C			
Temperature Stability	± 1.0°C			
Gross/ Net Weight (kg)	36/ 29 kg	48/ 42 kg	72/ 68 kg	79/ 77 kg
Power	500W	750W	1500W	2000W
Electrical Requirement	AC 220V 50/60Hz			