

Electronic Balance FA Series


- New Structure and Soft backlight LCD Display, giving rise to easy operation and outstanding performance
- Full automatic fault detection function ensures operation at safety
- Electromagnetic weighing sensor guarantees that products are weighed accurately
- Equipped with automatic calibration inside, which makes the operation of balance more simple and quick, the function more stable and the reading more accurate
- An external transformer is beneficial for accurate weighing and correct reading
- Internal RS 232 interface, enable to connect with computer, printer and other peripheral devices
- With such applications as counting, percentage, unit conversion and weighing, etc.



Specification of Electronic Balance FA Series

Model	ST-FA1004N	ST-FA1104N	ST-FA1604N	ST-FA2004N	ST-FA2004B	ST-FA2104N	ST-FA2104SN	ST-FA2204N	ST-FA3104N
Maximum Weighing (g)	100	110	160	200	200	210	210/60	220	310
Minimum Reading (mg)	0.1	0.1	0.1	0.1	0.1	0.1	1/0.1	0.1	0.1
Scale Size (mm)	ø 80								
Main Functions	Tare, calibration, auto-fault detection, overloading alarm, etc								
Output Interface	RS232C								
Overall Dimensions (mm)	356 x 330 x 215								
Power Sources	AC220V 50Hz								

Precision Balance



Specification of Precision Balance

Model	Capacity	Readability	Pan Size	Power Supply
ST-1002	100g	10mg	Ø 125mm	AC or DC
ST-2002	200g	10mg	Ø 125mm	AC or DC
ST-3002	300g	10mg	Ø 125mm	AC or DC
ST-4002	400g	10mg	Ø 125mm	AC or DC
ST-5002	500g	10mg	Ø 125mm	AC or DC
ST-10001	1000g	100mg	Ø 125mm	AC or DC
ST-10002	1000g	10mg	Ø 125mm	AC or DC
ST-12002	1200g	10mg	Ø 125mm	AC or DC
ST-20002	2000g	10mg	Ø 125mm	AC or DC
ST-1201	1200g	0.1g	Ø 125mm	AC or DC
ST-2001	2000g	0.1g	Ø 125mm	AC or DC
ST-3001	3000g	0.1g	145 x 178mm	AC or DC
ST-4001	4000g	0.1g	145 x 178mm	AC or DC
ST-5001	5000g	0.1g	145 x 178mm	AC or DC
ST-6001	6000g	0.1g	145 x 178mm	AC or DC
ST-15K	15kg	1g	290 x 280mm	AC or DC
ST-20K	20kg	1g	290 x 280mm	AC or DC
ST-30K	30kg	1g	290 x 280mm	AC or DC
ST-50K-1	50kg	1g	290 x 280mm	AC or DC
ST-50K	50kg	2g	290 x 280mm	AC or DC
ST-60K	60kg	2g	290 x 280mm	AC or DC