


Multi Position Stirring/ Stirring Hotplate



Specification of Multi Position Stirring/ Stirring Hotplate

Model	ST-AM6T	ST-AM6	ST-AM9T	ST-AM9	ST-AM12T	ST-AM12
Speed	0 ~ 1200 rpm					
Accuracy	5%					
Capacity	Capacity: up to 400ml each position					
Diameter of Function Rings	Diameter of function rings: 90mm each					
Number of Position	6		9		12	
Temperature Range	Ambient ~ 120°C					
Construction (Stirring only)	Stainless steel					
Construction (Stirring Hotplate)	Aluminium alloy					
With Hotplate	Yes	No	Yes	No	Yes	No
Heating Power	-	240W	-	360W	-	480W
Dimension (W x H x D) mm	278 x 55 x 298		278 x 55 x 388		278 x 55 x 478	

Stainless Steel Top Hotplate

- Temperature controller: thermostatic.
- Digital temperature display.
- Equipped with temperature probe.



Specification of Stainless Steel Top Hotplate

Model	ST-DB2A	ST-DB3A	ST-DB4A
Temperature Control Range	Ambient ~ 300°C		
Heating Power	500W	1000W	1500W
Stainless Steel Plate Size (mm)	150 x 200	200 x 300	300 x 400