

Walk-In Fume Hood

Model: ST-FHW-1500

- It's suitable for placing large instruments into the fume hood and operator is easy to enter into it.
- All steel covered with electrostatics spraying after folding, welding, polishing, acid cleaning, phosphorization, chemical resistance.
- Inner liner panel: imported from America, 6mm compact panel, acid and alkali resistance, temperature-resistance, anti-moisture.
- Guide plate: 6mm Compact panel, reasonable and knock-down guide structure
- Concentrated and effective discharge air, letting off different poisonous gas.
- Window: 6 mm thickness steel safe glass, balancing from elevation and subsidence, it can be stopped at any position by pulling.
- Lighting equipment: 30W antiseptis daylight lamp, lamp house is insulated with gas.
- Lightness on table not lower than 400 Lux.
- Operation panel: Lighting switch, ON/OFF switch two sets of 220V/10A versatile electric socket.
- External dimension: 1500 x 850 x 2300mm (W x D x H)

