

Rotary Evaporator ST-RE3000AA



- LCD touch screen controller for temperature setting, timer, water level alert display, low water cut-off, auto lift/down of glassware.
- Vertical glassware set consists of spiral condenser, vapour tube, stopcock for vacuum, delivery tube, PTFE sealing ring for maximum performance
- Temperature up to 90°C (Optional stainless steel bath with temperature up to 180°C without safety shield cover).

Specification of Rotary Evaporator ST-RE3000AA

Model	ST-RE3000AA
Rotation Speed	20~200 r.p.m.
Motor Rating	40 W
Heating Bath	Teflon bath 6L with protective safety shield cover
Heating Power	1.4 kW
Standard Evaporating Flask	1 L, joint 29/32 with clip
Standard Receiving Flask	1 L, joint 29/32 with clip
Temperature Range	Maximum 90°C (Optinal Stainless Steel bath with temperature up to 180°C without safety sheid cover)