




NUEVO

WIRE BASKETS

Made of AISI 304 stainless steel sheet with electro-bright finish. Designed to fit the Selecta autoclaves. (See page 54 to 67).

Part No.	Dimensions	Part No.	Dimensions
1004771	220 Ø x 140 mm high.	1000496	380 Ø x 280 mm high.
1000596	225 Ø x 240 mm high.	1000780	480 Ø x 220 mm high.
1000495	280 Ø x 200 mm high.		



ACCESSORY

Basket dividers, stainless steel AISI304, electrobright finish.

Part No.	Dimensions	For baskets, Part No.
1001217	4 compartment divider	1000495
1001218	6 compartment divider	1000496
1001222	7 compartment divider	1000780



WIRE BASKETS

Electropolished stainless steel AISI 304.

Part No.	Dimensions
1025115	210 Ø x 180 mm high
1025121	150 Ø x 120 mm high
1025124	240 Ø x 180 mm high

With cross dividers.

Part No.
1125115
1125121
1125124



STERILIZATION DRUMS

Made from AISI 304 stainless steel. Designed to fit the Selecta autoclaves. (See page 54 to 67.)

Part No.	Ø cm	Height cm
1002418	28	18
1002419	35	27
1002420	28	18

Doesn't match Selecta autoclaves.



STERILIZER CAN FOR PIPETTES


Made from Stainless steel AISI 304.

Part No.	height mm	Ø mm
1000829	400	72
1001126	400	130

STERILIZATION CONTAINER FOR PETRI DISHES

Made of AISI 304 stainless steel. Complete with Petri dish rack for up to 120 mm Ø.

Models	height / Ø mm	No. of plates
Part No.		
1001127	260 / 130	10
1001128	400 / 130	18



PETRI DISH RACK

Made of AISI 304 stainless steel, with electro-bright finish. For 10 dishes of 80 to 100 mm Ø.

Part No. 1000259



PIPETTE BASKET

Made of AISI 304 stainless steel with electro-bright finish. 90 mm Ø x 250 mm high.

Part No. 1000251



STERILIZATION CONTAINER FOR PIPETTES

Made from Stainless steel AISI 304. Lid interior and base lined with silicon to protect pipette end and tip.

Part No.	height mm	Ø mm
1001597	290	65
1001598	340	65