

Magnifying Lamp with Standard Table Bracket



Specification of Magnifying Lamp

Model	ST-MFL8606
Lens Size	5 inches
Diopter	5 diopter
Light Sources	22W circular fluorescent lamp

Portable Living Leaf Area Meter

- Measuring of leaves area without damaging the leaves.
- Easy to carry.
- Multiple measurements, such as leaf area, average area, length of leaf, width of leaf, thickness of leaf, perimeter and ratio of length to width.



Specification of Portable Living Leaf Area Meter

Model	ST-LLA1000
Test Unit	Milimeter, square milimeter
Precision	± or -2%
Length	≤ 1000mm
Width	≤ 160mm
Thickness	≤ 8mm
Stored Data	200 groups
Battery	6 V