

Portable pH Meter



- ST-PH4- Portable pH Meter with case and electrode

- ST-PH260- Portable pH Meter with case software and electrode

Specification of SASTEC Portable pH Meter			
Model	ST-PH260		ST-PH4
Parameters	pH/ mV (ORP)°C		pH/ mV (ORP)°C
Range	pH	0.00 ~ 14.00	
	mV	± 1800	± 1400
	Temperature	-5.0 ~ 105.0°C	-
Resolution	pH	0.01	
	mV	1	
	Temperature	0.1°C	-
Accuracy	pH	± 0.01	± 0.03
	mV	± 0.1%	± 0.2%
	Temperature	± 0.3°C	-
Temperature Compensation	Manual/ Auto: -5.0 ~ 105.0°C		Manual: 0.0 ~ 60.0°C
Calibration	May select 5 standard buffer solution for one and Two-point calibration		Two-point calibration
Data Storage	200 data with time-and-date		-
View, Delete and Print	Yes		-
Display	Backlit LCD		
Under-Voltage Indication	Yes		-
Output	RS-232		-
Software	1.1DC		-
Housing	IP65		-
Power	Four AA Batteries		Two AA Batteris
Case Outside Serial No.	WXS-A006-1		WXS-A004-1
Dimension (mm)	210 x 86 x 50		170 x 75 x 30
Weight (Kg)			0.5 Kg