

Class II Biosafety Cabinet

- Table top Model
- Automatic Sash
- Internal electrical socket
- Automatic UV cut off while sash is open
- 70% air circulation, 30% exhaust
- 304 Stainless Steel Interior
- HEPA Filter 99.999%
- Microprocessor controller



Specification of Class II Biosafety Cabinet

Model	ST-BSC700II
Overall Dimension (mm)	700 x 500 x 1200
Workzone Dimension (mm)	600 x 400 x 580
Noise Level	≤ 65dB(A)
Vibration	≤ 5um
Light	680 Lux
Power	800W
Power Supply	AC 220V 50Hz
Input Airflow	0.53 m/s
Output Airflow	0.33 m/s
Air Output	250m ³ /h

Eyewash Emergency Shower Station

Model: ST-WG7015F

This Floor Mounted Eyewash Emergency Shower combination of self draining shower head and eyewash with foot pedal.

- **Specification for Shower**
 Shower Head : Brass
 Stand Pipe : Stainless Steel (304 Grade)
 Activation : Stainless Steel pull rod operators a stay open stainless steel ball valves
- **Specification for Eyewash**
 Bowl : Stainless Steel (304 Grade)
 Eyewash Head : Rubber Cups
 Activation : Push flag operators a stay open stainless steel ball valves

