

Digital Ultrasonic Cleaner

- Double wall construction, made entirely of stainless steel.
- Draining stopcock for all models.
- 40kHz sound wave frequency, 59kHz available, factory preset.
- Digital LED display of timer, 0 to 80 minutes, temperature up to 90°C.
- Supply complete with basket and lid.



Specification of Digital Ultrasonic Cleaner

Model	ST-UB2200LT	ST-UB3200LT	ST-UB3300LT	ST-UB5200LT	ST-UB7200LT	ST-UB8200LT
Capacity	3L	4L	6L	10L	22L	22L
Frequency	40/59 kHz			28/40/59 kHz		
Power	120W	150W	180W	250W	500W	600W
Drainage	No		Yes		Yes	
Heating Function	Yes		Yes		Yes	
External Dimension (mm)	260 x 160 x 290	320 x 175 x 290	350 x 175 x 310	350 x 264 x 310	550 x 320 x 310	550 x 320 x 310
Internal Dimension (mm)	240 x 140 x 100	300 x 150 x 100	300 x 150 x 150	300 x 240 x 150	500 x 300 x 150	500 x 300 x 150

Dual Frequency Ultrasonic Cleaner

- Digital display for set and actual timer.
- Stainless steel interior and exterior.
- Equipped with lid, basket and drain.



Specification of Dual Frequency Ultrasonic Cleaner

Model	ST-300VDB
Exterior Dimension	320 x 264 x 320mm
Interior Dimension	300 x 240 x 150mm
Capacity	10L
Frequency	40/60 kHz
Power	300W
Adjustable Power	40 ~ 100%
Heating Power	400W
Temperature	10 ~ 80°C
Timer	1 ~ 480 min