

## Micropipettor, Adjustable Volume

- Reduce physical stress and discomfort, help avoiding repetitive strain injury.
- Application of modern innovative materials.
- Lightweight, ergonomic, low force design.
- Comply with ISO9001, 13485 and GLP guidelines.
- CE mark.
- Each pipettor has been factory-tested and certified at 22°C according to ISO8655/DIN12650.
- Each pipettors bears its own serial number and passes strict performance control attested by an individual QC certificate.
- One-handed operation for both right- and left-handers.
- The range of pipettors cover a volume range from 0.1ul to 5000ul.
- Easy to calibrate and maintain.
- Corrosion-resistant piston, tip cone and ejector collar.
- Three year warranty.



### Multi-channel pipettors

- Ejector works with dispensing head in any position
- 8- and 12-channel pipettors are available for standard 96-well plate.
- Dispensing head rotates for optimum pipetting convenience.
- Individual piston/tip cone assemblies allow easy repair and maintenance.
- Compound material tip cone design allows visual seal verification.

### Specification of Single Channel Pipette

Model	ST-770005	ST-770000	ST-770020	ST-770080	ST-770050	ST-770070	ST-770030	ST-770060	ST-770040	ST-770110
Volume (ul)	0.1 ~ 2.5	0.5 ~ 10	5 ~ 50	2 ~ 20	10 ~ 100	20 ~ 200	50 ~ 200	100 ~ 1000	200 ~ 1000	1000 ~ 5000
Increment (ul)	0.05	0.1	0.5	0.5	1	1	1	5	5	50

### Specification of Multi Channel Pipette

Model	8-Channel			12-Channel		
	ST-770210	ST-770220	ST-770240	ST-770310	ST-770320	ST-770340
Volume (ul)	0.5 ~ 10	5 ~ 50	50 ~ 300	0.5 ~ 10	5 ~ 50	50 ~ 300
Increment (ul)	0.1	0.5	5	0.1	0.5	5

## Micropipettor, Fixed Volume

- Reduce physical stress and discomfort, help avoiding repetitive strain injury.
- Application of modern innovative materials.
- Lightweight, ergonomic, low force design.
- Comply with ISO9001, 13485 and GLP guidelines.
- CE mark.
- Each pipettor has been factory-tested and certified at 22°C according to ISO8655/DIN12650.
- Each pipettors bears its own serial number and passes strict performance control attested by an individual QC certificate.
- One-handed operation for both right- and left-handers.
- Easy to calibrate and maintain.
- Corrosion-resistant piston, tip cone and ejector collar.
- Three year warranty.
- The range of fixed volume pipettes is designed for many different types of routine laboratory work.
- Using the same lightweight design and basic features of the Adjustable Volume range this product gives the user superior pipetting comfort and performance.
- Fixed volume pipette has white display to differentiate it from adjustable volume pipettes.



### Specification of Single Channel Pipette

Model	ST-722001	ST-722004	ST-722010	ST-722015	ST-722020	ST-722025	ST-722030	ST-722035	ST-722040	ST-722045	ST-722050	ST-722055
Volume (ul)	5	10	20	25	50	100	200	250	500	1000	2000	5000