



Fixed volume dispensing pipettes "Dispen-Fix"

WITH CALIBRATION CERTIFICATE WITH SERIAL NUMBER REFERENCE.

FEATURES

Fixed volume pipettes from 1 to 1000 µl, model dependent.

Simple and easy to use.

Discard tip ejection for all models.

Comes complete with calibration tool.

MODELS	Part No.	Capacity µl	Accuracy %	Precision %
Dispen-Fix	5904009	5	±2.0	< 0.8
Dispen-Fix	5904010	10	±1.6	< 0.5
Dispen-Fix	5904011	20	±1.0	< 0.3
Dispen-Fix	5904012	25	±1.0	< 0.3
Dispen-Fix	5904013	50	±0.8	< 0.3
Dispen-Fix	5904014	100	±0.8	< 0.25
Dispen-Fix	5904015	200	±0.8	< 0.15
Dispen-Fix	5904017	500	±0.8	< 0.2
Dispen-Fix	5904018	1000	±0.8	< 0.13



Original Comecta design



Variable volume dispensing pipettes "Dispen-Var"

WITH CALIBRATION CERTIFICATE WITH SERIAL NUMBER REFERENCE.

AUTOCLAVABLE UP TO 121 °C.

FEATURES

Variable volume from 1 to 5000 µl, model dependent.

Simple and easy to use adjustable dispense volume.

Digital set volume selection.


Discard tip ejection for all models up to 1000 µl.

Comes complete with calibration tool.

MODELS	Part No.	Capacity µl	Divisions µl	Precision %
Dispen-Var	5900010	1 - 10	0.1	< 1.5 - < 0.4
Dispen-Var	5900020	2 - 20	0.1	< 1.5 - < 0.8
Dispen-Var	5900100	10 - 100	1	< 1.0 - < 0.3
Dispen-Var	5900200	20 - 200	1	< 1.0 - < 0.5
Dispen-Var	5901000	100 - 1000	10	< 0.8 - < 0.28
Dispen-Var	5905006	1000 - 5000	10	< 0.4 - < 0.2

ACCESSORIES. PIPETTE TIPS. Accurate and precise, made of polypropylene.

Models N° image	Part No.	Range µl	Colour	Qty. tips x part No.	Suitable for pipette part No.	Volumes pipette µl
1	5490040	0.1- 10	Transparent	1000	5904009	5
					5904010	10
					5900010	1 - 10
2	5490041	2 - 20	Transparent	1000	5900020	2 - 20
					5904011	20
					5904012	25
3	5490042	5 - 200	Transparent	1000	5904013	50
					5904014	100
					5904015	200
4	5490044	101 - 1000	Blue	1000	5900100	20 - 100
					5900200	50 - 200
					5904017	500
5	5490047	1000 - 5000	Blue	500	5904018	1000
					5901000	100 - 1000
					5905006	1000 - 5000



Pipette rack.
Made of AISI 304 stainless steel. No. of places: 6 dosing pipettes up to 28 cm in length.
Part No. **1001215**