

Portable Bunsen Burner

- Precise and controlled
- Portable Instant Ignition
- Dual Use Locking Mechanism
- Refillable Butane Gas
- With dual flame, blue and yellow
- Flame length can be adjust
- Dual use locking system that will either lock the flame on for longer jobs or function as a safety lock.
- Gas inlet valve can be refill with butane gas.
- Ignition pad with instantly ignite.



Specification of Portable Bunsen Burner

Model	ST-K1000
Fuel Capacity	10 grams
Max Flame Temp (Blue)	1300°C
Max Flame temp (Yellow)	800°C
Continuos Burn Time Type	65 min (will vary with flame type & lenght)
Dimension (L x W x H) mm	143 x 67 x 68
Weight (grams)	348
Battery	1 AA Battery (Included)

IR Safety Gas Burner, Model ST-UF2000

- Ignite mode by touch-free IR sensing movement and a wireless foot pedal control.
- It has overheating and gas cut-off protection.
- Temperature at 1350°C on liquid.
- Portable, can be used with optional gas cartridge uses LPG (0.28~0.4mPa)
- Size (W x D x H): 100 x 130 x 95mm
- Weight: 2.5kg
- Can be connected to gas tank rechargeable battery

