



Upright freezers “Templow” S, M and L



APPLICATIONS

Preservation and storage of products that require to be preserved at low temperatures.

Upright freezer “Templow” S

**DIGITAL TEMPERATURE CONTROL,
TEMPERATURE RANGE FROM -18 °C TO -28 °C.
8 DRAWERS WITH INDIVIDUAL FREEZER COMPARTMENTS.**

FEATURES

Epoxy covered metallic external case.
Strong, molded, chemically resistant ABS interior.
Manual defrost with water evaporation system.
Door can be hung to open left or right.

CONTROL PANEL

Mains power ON/ OFF switch.
Temperature thermostat regulator.
Temperature indicator lamp.
Audible and visual alarms.
Push button alarm reset/stop.
Push button with quick freeze indicator.



Upright freezer “Templow” M

**TEMPERATURE RANGE FROM -9 °C TO -26 °C.
7 INDIVIDUAL FREEZER COMPARTMENTS WITH 14 BASKETS.**

FEATURES

Epoxy covered external case.
Strong, chemically resistant, ABS molded case.
Specially coated “Rilsan”, wire mesh baskets.
Manual defrost, with water evaporation system.
The door can be hung to open left or right.
Lockable door.

CONTROL PANEL

Temperature thermostat regulator.
Digital temperature display.



Upright freezer “Templow” L

**WITH FAN ASSISTED AIR CIRCULATION.
TEMPERATURE RANGE FROM -10 °C UP TO -36 °C.
COMES WITH 4 SHELVES.**

FEATURES

External and internal case made completely of AISI 310 stainless steel, with water sump and drain.
Chemically resistant.
Acclimatized hermetically sealed CFC free compressor.
Automatic defrost with electrical water evaporation system.
Regulable shelves.
Adjustable front level supports.
The door can be hung to open, left or right.
Door lock.

CONTROL PANEL

Mains power ON/ OFF switch.
Electronic temperature regulator.
Digital temperature display.



NEW DESIGN

MODELS	Part No.	Capacity litres	Range °C	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg	
TEMPLOW S	2101278	304	-18 to -28	150 44 42	184 60 63	110	74	Comes with 8 drawers
TEMPLOW M	2101279	513	-9 to -26	150 61 57	186 75 76	345	111	*Comes with 14 baskets
TEMPLOW L	2101280	600	-10 to -36	155 51 65	215 70 83	600	160	*Comes with 4 shelves

Comes with steel PVC laminated shelves and baskets.