

Table Top Freeze Dryer



- SASTEC Freeze Dryer uses environmental friendly CFC Free Refrigerants.
- SASTEC Freezer Dryer uses the force air cooling condensation type refrigerationsystem.
- Transparent drying chamber, it has four Stainless Steel shelves to accept 22mm vials. Drying chamber has 8 ports for freezing flasks.
- Four shelves can accomodate
- 135 pieces of 22mm vials/ tubes.
- Wide mouth cold trap can have pre-freeze function.
- Low noise compressor.
- All Stainless Steel construction against corrosion.
- Compact design with stainless steel trolley.
- LCD for display of setting and operating parameters.
- CFC-Free refrigeration system to cool collector to 50°C (Optional for -80°C)

Specification of Table Top Freeze Dryer

Model	ST-FD-10B-MPT
Vacuum	< 15 Pa
Water Removal	3 L/25 Hours
Dimension (L x W x H) mm	580 X 560 X 340
Gross Weight	85 kg
Holding Ice Capacity	Up to 1.2 L