

Medical/Vaccine Refrigerator



ST-MD 120
 ST-MD 170
 ST-MD 210



ST-MD 280
 ST-MD 300
 ST-MD 340



ST-MD 400
 ST-MD 560
 ST-MD 600



ST-MD 800



ST-MD 1200

- The Medicine/Vaccine Refrigerator is the specific medicine cold-storage equipment. It can be also used for storing reagent, vaccine, biological product, etc. It is especially for hospital, pharmacy, pharmaceutical factory, epidemic station, blood bank station, health center and laboratory.
- **Precise temperature control system**
 The high-accuracy microcomputer-controlled temperature system maintains the temperature constancy at 2~8°C precisely. The temperature set-value can be adjusted from 0°C~14°C, and the temperature display precision is 0.1°C, the resolution is 0.1°C and the adjusting increment is 0.1°C. The value is to be displayed through three-digit high-brightness figure. The system has special high-resolution temperature calibrating function. That means the set-value figure can be displayed and calibrated at 0~14°C. The calibrating increment is 0.1°C. The forced cold-air circulation system ensure the temperature uniformity inside the chamber.