



## Ultrasonic cleaning baths "Ultrasons" without heating



### CONTROL PANEL

Timer switch for 0 -15 minutes or continuous operation. Indicator lamp when the unit is in operation.

### MODELS

Part No.	Capacity litres	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Generator power W	Timer 0-15'	Waste drain	Weight Kg
3000512 1	1	6 14 14	20 20 18	50	Yes	No	2,5
3000683 2	2,6	9 23 13	22 29 19	100	Yes	No	4
3001208 8	4	10 30 14	24 33 18	100	Yes	Yes	7
3000513 3	6	15 30 14	30 36 20	150	Yes	Yes	9,5
3000837 5	6,5	14 24 Ø	30 27,5 Ø	150	Yes	No	8
3000617 4	9	15 50 14	30 56 20	200	Yes	Yes	10,7
3000514 6	9	15 30 23	30 36 30	200	Yes	Yes	13
3000515 7	20	15 50 30	30 56 36	400	Yes	Yes	22

### COMPLEMENT

**Detergent "BIO-SEL"**  
Bio-degradable, phosphate free. Excellent cleaning power just by mixing only 20ml per litre of water. BIO-SEL, in a 2% solution has a 9.7 pH level. Low toxicity. This detergent disperses, emulsifies, segregates, dissolves, suspends and breaks down dirtiness.

**1000025** Capacity 1 litre.  
**1000026** Capacity 4 litres.  
**1000852** Box of 12 x 1 litre bottles.  
**1000853** Box of 4 x 4 litres bottles.



### ACCESSORIES

**Baskets and Lids.**  
Made of AISI 310 stainless steel.



Cleaner model Part No.	Part No	BASKETS Height / Width / Depth (cm)	LID WITH CONCENTRIC REDUCTION RING Part No.	Number of positions	LID WITH HANDLE Part No.
3000512	6005123	5 11 11	6005122	1	6005121
3000683	6006833	7 20 12	6006832	1	6006831
3001208	6001209	7 26 13	-	-	6005131
3000513	6005133	12 26 13	6005132	2	6005131
3000837	6008373	13 23 Ø	-	-	6008371
3000617	6005144	12 46 12	6005134	3	6005124
3000514	6005143	13 29 23	6005142	2	6005141
3000515	6005153	13 49 29	6005152	6	6005151



## Ultrasonic Bath "Ultrasons-Digit"

5 CYCLES ELECTRONIC DIGITAL TIME CONTROL  
CAPACITY 1.4 LITRES

### FEATURES

Internal tank made from stainless steel, with minimum and maximum level indicator.  
External case, basket and lid made from PVC.  
Maximum cleaning power, frequency 42 KHz

### CONTROL PANEL

Pushbutton start/ stop  
Timer 5 programs selection in seconds:  
90" - 180" - 280" - 380" and 480".  
Wide application range from jewellery, watches, optical components. Dentistry and general laboratory use.

### MODEL

Part No.	Capacity litres	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
5900650	1.4	6 13 16	15 17 22	70	1.75

Comes complete with basket and see through transparent PVC lid.



OPTIC  
**ivymen**  
SYSTEM

Model "Elegant".