

Vortex Mixer ST-MXS



- Orbital shaking movement with a 4 mm diameter mixes samples quickly and easily with continuous or touch operation.
- Sturdy steel base with silicon feet to absorb vibration at high speeds.
- Stepless speed regulation of motor. Eccentric with oil-less ball bearings.
- Can accommodate Eppendorf tubes, microtiter plates and Erlenmeyer flask when used with optional adapters and inserts.
- Dimension: 127 x 130 x 160mm
- Weight: 3.5kg
- Adjustable speed 0~2500 rpm
- Warranty: 1 year

Parts & Accessories of Vortex Mixer ST-MXS

Model	Description
ST-MXS-700130	Universal adapter for various rubber foam inserts
ST-MXS-700131	Test tube insert for 48 test tubes, 6 mm diam.
ST-MXS-700132	Test tube insert for 18 test tubes, 10 mm diam.
ST-MXS-700133	Test tube insert for 12 test tubes, 12 mm diam.
ST-MXS-700134	Test tube insert for 8 test tubes, 16 mm diam.
ST-MXS-700135	Test tube insert for 8 test tubes, 20 mm diam.