

ACCESSORIES for the TECTRON-BIO, TECTRON-200, DIGITERM-100 and, DIGITERM-200. Made of AISI 304 stainless steel.



1. Gable lid.

Part No. **3001295**

2. Lifting rack support. Capacity 3 tube racks.

Part No. **1001296**

3. Tube racks for lifting rack support.

Part No.	For Ø tubes mm	Capacity Rack tubes	Height/Width/Depth (exterior) cm
1001202	13	36	8 8.2 23.6
1001203	16	24	8 8.2 23.6
1001204	20	24	8 8.2 23.6



Immersion thermostat for baths "Digit-Cool"

DIGITAL SELECTOR AND TEMPERATURE CONTROL.

ADJUSTABLE TEMPERATURE FROM 0 °C TO 100 °C.

STABILITY: 100 °C ±0.1 °C, AND AT 0 °C ±0.05 °C. SET ERROR: ±0.1 °C. RESOLUTION: 0.1 °C.



SAFETY:

STANDARD DIN 12879.2

SAFETY OVER TEMPERATURE THERMOSTAT WITH ADJUSTABLE MANUAL RESET LOW LIQUID LEVEL SENSOR PROTECTION.

Autonomous unit with cooling system

Compressor unit built-in

Features and control panel see page 98.

FEATURES

Maximum volume to achieve the working temperature: 20 litres.

Minimum depth of the elements with contact with the liquid: 15 cm.

Hermetically sealed compressor.

RS-232 out put to a computer control or printer.

Portable control and stirring unit with lifting handle.



MODEL

Part No.	Control Range °C	Height / Width / Depth (exterior) cm	Pump Pressure / Flow H ₂ O	Power W	Power H.P.	Weight Kg
3001373	0 +100	45 21 51	150 mbar / 12 l/min	1460	1 / 5	22

ACCESSORIES

Stainless steel baths.

Double skinned with an AISI 304 stainless steel exterior with a bath of AISI 316 stainless steel interior.

Supplied with a drain tap and lifting handles.



1. Lifting rack support. .

Capacity 3 tube racks. Part No. **1001296**

2. Tube racks for lifting rack support.

Part No.	For Ø tubes mm	Capacity Rack tubes	Height / Width / Depth (exterior) cm
1001202	13	36	8 8.2 23.6
1001203	16	24	8 8.2 23.6
1001204	20	24	8 8.2 23.6

MODELS

Part No.	Capacity litres	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Weight Kg
6000392	20	15 48 30	20 55 35	6.6
6000393	27	20 48 30	25 55 35	7.6



Ring set lid 105, 80, 60 and 37 mm Ø reduction rings and a check thermometer location hole in the lid.

Capacity 4 places. Part No. **1001374**

Digitcool thermostat with bath and 4 place ring set lid.