



## Centrifuge “Centro-4-BL”

WITH SWINGOUT HEAD FOR 4 TUBES OF 75 X 13 mm VAC.  
WITH FREE FROM MAINTENANCE INDUCTION DRIVE MOTOR.



Control panel features (see page 101).

Centrifuge kit for separating plasma rich platelets from whole blood.



ABS plastic moulded exterior.



### MODEL

Part No.	Height / Width / Depth (exterior) cm	Power W	Max. R.p.m.	Max. R.c.f. (xg)	Weight Kg
<b>7002560</b>	<b>20 28 34</b>	<b>170</b>	<b>5000</b>	<b>2907</b>	<b>6</b>

Rotor comes without tubes.



## Centrifuge “Centro-8-BL”

WITH ANGLE ROTOR FOR 8 X 15 ml TUBES, OR 7 ml TUBES OR TUBES SIZES 75 X 13 mm VAC (WITH ADAPTERS).  
WITH FREE FROM MAINTENANCE INDUCTION DRIVE MOTOR.

Control panel features (see page 101).



ABS plastic moulded exterior.



### MODEL

Part No.	Max. Volume ml.	Tube capacity	Height / Width / Depth (exterior) cm	Power W	Max. R.p.m.	Max. R.c.f. (xg)	Weight Kg
<b>7002356</b>	<b>120</b>	<b>8 x 15 ml</b>	<b>20 28 34</b>	<b>170</b>	<b>5000</b>	<b>2515</b>	<b>6.5</b>

Comes complete with 8 x 15 ml glass tubes.

### ACCESSORIES

**Bucket adapter for 7 ml tubes.**  
Rotor capacity 8 adapters.  
Part No. **7000700**

**Bucket adapter for 75 x 13 mm VAC. tubes.**  
Rotor capacity 8 adapters.  
Part No. **1001557**

**SPARES. Tubes** (see page 118).