



Centrifuge “Centromix II-BL”

WITH A CHOICE OF ROTORS AND INTERCHANGEABLE ADAPTERS.
WITH FREE FROM MAINTENANCE INDUCTION DRIVE MOTOR.



NEW

Control panel features (see page 101).



External DUR-AL case, with AISI 304 stainless steel bowl.




MODEL

Part No.	Max. Volume ml	Max. Number of tubes	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7002555	120	4 x 25 ml	23 30 35	180	12

ACCESSORIES

Rotors and interchangeable buckets and adapters.

Part No. head	Rotor type	Part No. bucket	Total max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
7000550	Angle	-	120 ml	8 x 15 ml	16.5 x 100	4400	2063
7001065	Swing out	-	-	-	-	-	-
	Swing out	7001066	100 ml	1 x 25 ml	24 x 100	4400	2601
	Swing out	7001067	60 ml	1 x 15 ml	16.5 x 100	4400	2601
7001065	Swing out	7001068	28 ml	1 x 7 ml	12 x 100	4400	2601
	Swing out	7001573	40 ml	2 x 5 ml VAC	13 x 75	4400	2163
	Swing out	7001070	60 ml	VAC	16 x 100	4400	2601
	Swing out	7001072	-	undrilled	-	4400	2601



7000550

7001065



7001066

7001067

7001068

7001573

7001070

7001072

Nota: All buckets and adapters come without tubes unless stated, see page 118 for tubes and accessories.



Centrifuge “Mixtasel-BL”

MAINTENANCE FREE INDUCTION DRIVE MOTOR.



NEW

Control panel features (see page 101).



External DUR-AL case, with AISI 304 stainless steel bowl.



Polyvalent

With rotors and interchangeable adapters

MODEL

Part No.	Max. Volume ml	Max. Number of tubes	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7002575	400	4 x 100 ml	29 39 44	210	22