



Centrifuge "Centromix II-BL"

WITH A CHOICE OF ROTORS AND INTERCHANGEABLE ADAPTERS.
FREE MAINTENANCE INDUCTION DRIVE MOTOR.

This model can be marked as
medical device



Control panel features (see page 111).



MODEL

Part No.	Max. Volume ml	Max. Number of tubes	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7002555	120	4 x 25 ml	23 30 35	180	12

External DUR-AL case, with AISI 304 stainless steel bowl.

ACCESSORIES

Rotors and interchangeable buckets and adapters.

Part No. head	Rotor type	Part No. bucket	Total max. Volume (4 adapters)	Bucket/Adapter capacity	Ø x Length mm	Max. R.p.m.	Max. R.c.f. (xg)
7000550	Angle	-	120 ml	8 x 15 ml	16.5 x 100	4400	2063
7001065	Swing out	-	-	-	-	-	-
	Swing out	7001066	100 ml	1 x 25 ml	24 x 100	4400	2601
	Swing out	7001067	60 ml	1 x 15 ml	16.5 x 100	4400	2601
7001065	Swing out	7001068	28 ml	1 x 7 ml	12 x 100	4400	2601
	Swing out	7001573	40 ml	2 x 5 ml VAC	13 x 75	4400	2163
	Swing out	7001070	60 ml	VAC	16 x 100	4400	2601
	Swing out	7001072	-	undrilled	-	4400	2601



Note: All buckets and adapters come without tubes unless stated, see page 127 for tubes and accessories.



Centrifuge "Mixtasel-BL"

MAINTENANCE FREE INDUCTION DRIVE MOTOR.



Polyvalent with rotors and interchangeable adapters

Control panel features (see page 111).



NEW DESIGN

External case of DUR-AL alloy with ABS plastic coated front AISI 304 stainless steel bowl.

MODEL

Part No.	Max. Volume ml	Max. Number of tubes	Height / Width / Depth (exterior) cm	Power W	Weight Kg
7002575	400	4 x 100 ml	29 39 44	210	22