

Upright Two Door Freezer

- Temperature Range: -10 to -40°C
- Microprocessor temperature controller, LCD display
- Audio/Visual alarm for High/Low temperature
- 90mm CFC-free insulation
- 2 units of compressor (Danfoss) and EBM motor fan.
- Double sealing with lock
- Painted metal external and ABS material internal.
- Three shelves each for upper and lower chamber.
- Mounted on castors.
- CFC-free refrigerant R507.
- Certification: ISO 9001, ISO 14001, ISO 13485
- Power supply: 220V 50Hz (110V 60Hz is available)
- Warranty period: 2 years on all parts.



ST-FF40-208

Specification of Upright Two Door Freezer		
Model	ST-FF40-208	
Type	Upper	Lower
Capacity (L)	44	109
External Dimension (W x D x H) mm	623 x 604 x 1775	
Internal Dimension (W x D x H) mm	397 x 368 x 781	397 x 668 x 414
Input Power (W)	370	
Weight (kg)	95	