


# Organomation™

Laboratory Evaporators and Extractors



**Blowdown Nitrogen  
Evaporators for  
Sample Vials  
Large Volume  
Evaporators with  
Solvent Collection**



**Your Partners  
in Sample  
Preparation  
since 1959**



**Organomation  
Associates, Inc.**  
266 River Road West  
Berlin, MA 01503-1699  
USA

978.838.7300 • [organomation.com](http://organomation.com)

# Solvent Extraction

You are in the lab and need to run multiple extractions at once. Between the individual Soxhlet extractors, heating mantels and tubing, you will spend more time assembling and monitoring your apparatuses than actually performing chemistry.

The ROT-X-TRACT product line eliminates chaos and wasted time by allowing for up to 8 or 10 extractions to be performed simultaneously in a single bath with a pair of water manifolds.

ROT-X-TRACT-S



Rotary Solid-Liquid  
Soxhlet Extractor

ROT-X-TRACT-L



Rotary Liquid-Liquid  
Extractor

ROT-X-TRACT-LC



Corning Accelerated  
One-Step  
Liquid-Liquid Extractor

## Advantages:

- **Conserve valuable bench space:** All samples are arranged in a circle.
- **Easy sample access:** Each instrument rotates allowing all samples to be accessed from the front.
- **Fewer connections:** One water supply line in, one drain line out.